

Recommended Cap Sizes per Age, updated reference table

Cap size	Age	Can comfortably accommodate this number of electrodes
20	Preterm week 22	19
22	Preterm week 24	19
24	Preterm week 26	19
26	Preterm week 28	19
28	Preterm week 30	19
30	Preterm week 32	32
32	Preterm week 34	32
34	Preterm week 37	32
36	Newborn	32
38	1 month	32
40	3 months	32
42	5 months	64
44	7 months	64
46	1 year	64
48	2 years	64
50	3-4 years	64
52	5-10 years	128
54	11-14 years, Adult small female	128
56	15-17 years, Adult medium (female)	128
58	18 years, Adult medium (male)	128
60	Adult large male	128
62	Adult very large	128
64	Adult very very large	128