

**Press Release**

March 2011, Volume 38

Dear Customers, Dear Friends of Brain Products,



we at Brain Products care greatly about even the tiniest details of our products. We truly believe that it is not sufficient merely to provide you with a functioning system. In our view, our products must be both good AND user-friendly. That is why I am very happy to inform you that a new version of the actiCAP ControlBox is now available.

The new design of the ControlBox II makes it extremely easy to replace or add electrodes, which was a tedious chore with the ControlBox I.

Electrodes are no longer permanently attached to the ControlBox as in the previous version. All you have to do now to add or replace electrodes is to plug them into the connectors located on the front

panel of the actiCAP ControlBox II. After rebooting, the number of electrode sets connected to the actiCAP system will automatically be detected by the ControlBox.

Using the ControlBox II requires actiCAP ControlSoftware 1.2.2.0, which can be downloaded gratis from our website (<http://www.brainproducts.com/downloads.php?kid=8&tab=4>). The previous version of the ControlBox is also compatible with actiCAP ControlSoftware 1.2.2.0.

Customers ordering on or after March 15, 2011 will automatically be sent the new actiCAP ControlBox.

Naturally, existing actiCAP customers can upgrade their ControlBox I to a ControlBox II for a very modest exchange fee. Please contact either us or your local dealer for more information.

Please note that you will need to send back the complete actiCAP system to upgrade your ControlBoxes.

We hope you will enjoy reading this issue of our Press Release!

Pierluigi Castellone  
Brain Products' CEO

**Contents of this issue**

**IN THE FOCUS**  
Brain Products is moving forward in the BCI field 2

**User Research**  
Multi-stability and non-linearity of large-scale cortical rhythms 2

**Did you know...!?**  
Towards passive BCI (3) 6

**Support Tip**  
Introducing the new „All-In-One-Installer“ for BrainVision Analyzer 2 8

**News in brief**  
Workshops 8

**Products in Practice**  
Exploring the EEG correlates of meditation 9

**Products in Practice**  
EEG in weightless conditions 10

**Inside Regulatory Affairs**  
Brain Products achieves certification for Product and Service Quality 11

**Young Scientist Award**  
UCL-CNT Young Investigator Award in Neuroimaging Techniques 12

**Customer Survey**  
Brain Products Online Survey 12

**Brain Products Projects**  
KALOMED 13

**Brain Products Inside**  
Who is who - Stefanie Rudrich 13

**Brain Products Inside**  
Who is who - Birgit Trogisch 13


**News in brief**  
Downloads, Programs and Updates 14

**News in brief**  
Conferences 14

Contact & Imprint 14



ControlBox Version I



ControlBox Version II, ribbon cable for the connection to the amplifier (connectors on the lower part of the ControlBox) and splitter box with actiCAP electrodes (connectors on the upper part of the ControlBox).